

Specifications

# Shimmer3R ECG/EMG Unit



User Guides  
Sample Data  
Case Studies  
Software  
Expansion Boards

## Key Features

5-wire, 4-channel ECG, measuring bipolar limb leads and user's choice of V1 - V6  
Measure 2 channels of EMG data with a common reference electrode  
Software configurable right-leg drive for common-mode interference rejection  
Software configurable amplifier gain  
Software configurable data rate  
Respiration demodulation capability on-chip  
Lead-off detection capability on-chip  
Test signal on-chip for validation purposes  
EEPROM storage device (on the ExG daughterboard) enables expansion board detection and identification, as well as 2032 bytes of data storage available to user

## Introduction

The ECG/EMG unit, previously also known as the ExG unit, provides a configurable digital front-end, optimised for the measurement of physiological signals, for example 5-wire ECG (Electrocardiography) and 2-channel EMG (Electromyography). Compatible with the Shimmer3 platform, the ECG/EMG unit also boasts best data quality with integrated 10 DoF inertial sensing via accelerometer, gyroscope, magnetometer and altimeter, each with selectable range. The ECG/EMG unit also provides highly accurate and scientifically reliable raw data to allow complete control over capture and interpretation of sensed data in real-time.

## Product Overview

While addressing the challenges of mobility, the ECG unit can record the pathway of electrical impulses through the heart muscle, and can be recorded on resting and ambulatory subjects.

When configured for EMG, the unit can also measure and record the electrical activity associated with muscle contractions, assess nerve conduction, muscle response in injured tissue, activation level, or can be used to analyse and measure the biomechanics of human movement.

## Applications

Atrial fibrillation  
Premature ventricular contraction  
Heart function monitoring  
Abnormal rhythm detection and alert  
Biomechanics, muscle activity, gait and posture disturbance  
Fatigue analysis  
Sports technique, performance and medicine  
Neuro Rehabilitation  
Tremor Analysis  
Veterinary Science  
Orthopaedics  
Thoracic Bioimpedance (Respiration)

## Technical Specifications

|                                  |                                                                                                                                                                                                         |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gain                             | Software configurable (1, 2, 3, 4, 6, 8, 12)                                                                                                                                                            |
| Data Rate                        | Software configurable<br>(125, 250, 500, 1000, 2000, 4000, 8000 SPS)                                                                                                                                    |
| Input Differential Dynamic Range | Approx 800 mV (for gain = 6)                                                                                                                                                                            |
| Bandwidth                        | 8.4 kHz                                                                                                                                                                                                 |
| Ground                           | Wilson Type Driven Ground                                                                                                                                                                               |
| Input Protections                | ESD and RF/EMI filtering; Current limiting; inputs include defibrillation protection (survive only, not repeat)                                                                                         |
| Connections                      | EMG: Input Ch1N, Input Ch1P, Input Ch2N, Input Ch2P, Reference (Ref)<br>ECG: Input RA, Input LA, Input LL, Input Vx, Reference (RL)<br>All Hospital-Grade 1mm Touchproof IEC/EN 60601-1 DIN42-802 jacks |
| Weight                           | 31 grams                                                                                                                                                                                                |
| Dimensions                       | 65 x 32 x 12 mm                                                                                                                                                                                         |
| EEPROM Memory                    | 2048 bytes                                                                                                                                                                                              |

## Shimmer3R Specifications

|                          |                                                  |
|--------------------------|--------------------------------------------------|
| Processing               | STM32U5A5 Arm Cortex-M33 (160 MHz, 16-bit)       |
| Communication            | Classic Bluetooth and BLE (Vela IF820)           |
| Storage                  | Integrated MicroSD card slot (8GB card supplied) |
| Battery                  | 450mAh rechargeable Li-ion                       |
| Integrated 3-Axis Accel. | LIS2DW12 and the LSM6DSV                         |
| Accel. Range             | ±2g , ±4g , ±8g , ±16g                           |

|                   |                                                                                                                                                                                                                                                |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ECG Module</b> | <p>5-wire, 4- channel ECG, measuring bipolar limb leads and user's choice of V1 - V6</p> <p>Digital interface includes test signal generation for validation purposes</p> <p>Respiration demodulation from ECG data and lead-off detection</p> |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|                   |                                                                                                                  |
|-------------------|------------------------------------------------------------------------------------------------------------------|
| <b>EMG Module</b> | <p>Two channels of EMG data</p> <p>Digital interface includes test signal generation for validation purposes</p> |
|-------------------|------------------------------------------------------------------------------------------------------------------|

|                            |                                                                                                                                                                                                                           |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Supporting Software</b> | <p>Shimmer ConsensysPRO &amp; ConsensysBASIC</p> <p>Shimmer LabVIEW Instrument Driver</p> <p>Shimmer MATLAB Instrument Driver</p> <p>Shimmer Java / Android API</p> <p>Shimmer C# API</p> <p>Shimmer 9DoF Calibration</p> |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|                                              |                                                                                                                                 |
|----------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Supporting Hardware &amp; Accessories</b> | <p>Biophysical Leads</p> <p>Straps, Documents, Charging Dock/Base, Case</p> <p>Electrodes</p> <p>Shimmer3 Calibration Stand</p> |
|----------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|